



**memmert**  
Experts in Thermostatics



## ADDITIONAL ACCESSORIES

100% ATMOSAFE. MADE IN GERMANY.

[www.memmert.com](http://www.memmert.com) | [www.atmosafe.net](http://www.atmosafe.net)

## Manufacturer and customer service

Memmert GmbH + Co. KG  
Willi-Memmert-Straße 90–96  
D-91186 Büchenbach  
Germany

Phone: +49 (0)9122 925-0  
Fax: +49 (0)9122 14585  
E-mail: [sales@memmert.com](mailto:sales@memmert.com)  
Internet: [www.memmert.com](http://www.memmert.com)

### Customer service:

Service hotline: +49 (0)9171 9792 911  
Service fax: +49 (0)9171 9792 979  
E-mail: [service@memmert.com](mailto:service@memmert.com)

### Shipping address for repairs:

Memmert GmbH + Co. KG  
Customer service  
Willi-Memmert-Str. 90-96  
DE-91186 Büchenbach  
Germany

Please contact our customer service before sending appliances for repair or before returning equipment, otherwise, we have to refuse acceptance of the shipment.

© 2016 MEMMERT GmbH + Co. KG

D24124 | Date 11/2016

We reserve the right to make changes

## About this manual

---

### Purpose and target audience

This manual describes the optional additional accessories for Memmert Generation 2012 appliances. It is intended for use by trained personnel of the owner, who have the task of operating and/or maintaining the respective appliance.





### Other documents that have to be observed:

- ▶ the operating manual of the respective appliance
- ▶ For operation of the appliance with MEMMERT AtmoCONTROL, observe the separate software manual





## Contents

About this manual	3
Contents	4
Floating switching contact ALARM (opt. H6)	5
Floating switching contact setpoint reached (SP) (opt. H5)	6
Freely programmable switching contact (A – D)	7
4–20-mA current loop interface for temperature	8
4–20-mA current loop interface for humidity	8
Electric door locking mechanism	9
Freely positionable Pt100 temperature sensor (option H8)	9
HPP constant climate chamber with compressed air dehumidification	10
Gas flushing	11







## Floating switching contact ALARM (opt. H6)

Relay/ LED	Plug assignment	Functional description	Miscellaneous
Relay off  Test LED red  	Alarm 	Contact 2-3 is closed in case of the following errors: <ul style="list-style-type: none"> <li>▶ Loss of voltage</li> <li>▶ Overtemperature</li> <li>▶ Undertemperature (plus controller only)</li> <li>▶ Humidity alarm</li> <li>▶ Mechanical Temperature limiter TB</li> <li>▶ Fan speed alarm (only for optional fan speed monitoring)</li> <li>▶ Error of sensor PT100</li> </ul>	Switching capacity max. 2 A max. 24 Volt
Relay on  Test LED green  	Alarm 	Contact 1-2 is closed <ul style="list-style-type: none"> <li>▶ Appliance switched on and in order</li> </ul>	Switching capacity: max. 2 A max. 24 V


## Floating switching contact setpoint reached (SP) (opt. H5)

Relay/ LED	Plug assignment	Functional description	Miscellaneous						
Relay off  Test LED red  	SP 	Contact 2-3 is closed  ▶ Setpoint not reached  ▶ The actual temperature value is outside of a set tolerance band around the defined setpoint  <table><tr><th>Appliance type</th><th>Tolerance band</th></tr><tr><td>IN, IF, IPP, HPP, INCO2, HCP, ICP, ICH, IPS</td><td>dT ≥ 0.5 K</td></tr><tr><td>UN, UF, SN, SF, CTC, TTC</td><td>dT ≥ 2.0 K</td></tr></table>	Appliance type	Tolerance band	IN, IF, IPP, HPP, INCO2, HCP, ICP, ICH, IPS	dT ≥ 0.5 K	UN, UF, SN, SF, CTC, TTC	dT ≥ 2.0 K	Switching capacity max. 2 A max. 24 Volt
Appliance type	Tolerance band								
IN, IF, IPP, HPP, INCO2, HCP, ICP, ICH, IPS	dT ≥ 0.5 K								
UN, UF, SN, SF, CTC, TTC	dT ≥ 2.0 K								
Relay on  Test LED green  	SP 	Contact 1-2 is closed  ▶ Setpoint reached  ▶ The actual temperature value is within a set tolerance band around the defined setpoint  <table><tr><th>Appliance type</th><th>Tolerance band</th></tr><tr><td>IN, IF, IPP, HPP, INCO2, HCP, ICP, ICH, IPS</td><td>dT &lt; 0.5 K</td></tr><tr><td>UN, UF, SN, SF, CTC, TTC</td><td>dT &lt; 2.0 K</td></tr></table>	Appliance type	Tolerance band	IN, IF, IPP, HPP, INCO2, HCP, ICP, ICH, IPS	dT < 0.5 K	UN, UF, SN, SF, CTC, TTC	dT < 2.0 K	Switching capacity: max. 2 A max. 24 V
Appliance type	Tolerance band								
IN, IF, IPP, HPP, INCO2, HCP, ICP, ICH, IPS	dT < 0.5 K								
UN, UF, SN, SF, CTC, TTC	dT < 2.0 K								

## Freely programmable switching contact (A – D)

Relay/ LED	Plug assignment	Functional description	Miscellaneous
Relay off  Test LED red  	Out A (B,C,D) 	Programming of the freely programmable switching contacts is done using the AtmoCONTROL software.  Up to 4 (for single-phase appliances max. 2) floating switching contacts can be switched programme dependently. For switch setting "open" ► Contact 1-2 open ► Contact 2-3 closed	Switching capacity max. 2 A max. 24 Volt  
Relay on  Test LED green  	Out A (B,C,D) 	For switch setting "close" ► Contact 1-2 closed ► Contact 2-3 open	Switching capacity: max. 2 A max. 24 V  


## 4–20-mA current loop interface for temperature

Plug assignment	Appliance	Range	4 mA	12 mA	20 mA
4-20mA / °C 	IN / IF	0 ... + 90 °C	0 °C	45 °C	90 °C
	IN <sup>PLUS</sup> /IF <sup>PLUS</sup> (with Steri function)	0 ... + 90 °C	0 °C	45 °C	90 °C
	SN / SF	0 ... + 260 °C	0 °C	130 °C	260 °C
	UN / UF	0 ... + 310 °C	0 °C	155 °C	310 °C
	IPP / HPP	–10 ... + 80 °C	– 10 °C	35 °C	80 °C
	ICP / ICH	–20 ... + 70 °C	– 20 °C	25 °C	70 °C
	CTC	–50 ... + 200 °C	– 50 °C	75 °C	200 °C
	Optional	0 ... + 70 °C	0 °C	35 °C	70 °C
	Optional	0 ... + 80 °C	0 °C	40 °C	80 °C
	Optional	0 ... + 100 °C	0 °C	50 °C	100 °C
	Optional	0 ... + 300 °C	0 °C	150 °C	300 °C
	Optional	20 ... + 90 °C	20 °C	55 °C	90 °C
	Optional	20 ... + 100 °C	20 °C	60 °C	100 °C
	Optional	20 ... + 200 °C	20 °C	110 °C	200 °C
	Optional	20 ... + 260 °C	20 °C	140 °C	260 °C
	Optional	20 ... + 300 °C	20 °C	160 °C	300 °C
	Optional	20 ... + 310 °C	20 °C	165 °C	310 °C

R resistance: max.  $2,5 \text{ V} @ 20 \text{ mA} = 125 \text{ Ohm}$

For errors 0 mA output.

## 4–20-mA current loop interface for humidity

Plug assignment	Functional description	Miscellaneous								
<div>4-20mA / °C</div> 	Current loop interface 4-20 mA	R resistance:								
	<table><tr><td>Range</td><td>4 mA</td><td>12 mA</td><td>20 mA</td></tr><tr><td>0...100 % rh</td><td>0 % rh</td><td>50 % rh</td><td>100 % rh</td></tr></table>	Range	4 mA	12 mA	20 mA	0...100 % rh	0 % rh	50 % rh	100 % rh	Max. 2.5V@ 20mA = 125 Ohm
	Range	4 mA	12 mA	20 mA						
0...100 % rh	0 % rh	50 % rh	100 % rh							




# Electric door locking mechanism


## Functional description

Programming of the electric door locking mechanism is done using the AtmoCONTROL software. You can set the door locking mechanism to be electrically locked and unlocked at any stage of the programme.

Setting "close":  
Door locking mechanism  
electrically locked




Setting "open":  
Door locking mechanism  
electrically unlocked





## Freely positionable Pt100 temperature sensor (option H8)

The Pt100 temperature sensor can be flexibly positioned in the interior or in the chamber load to measure temperatures locally (a maximum of 3 additional sensors is possible). The individual temperatures measured are logged in the integrated data logger and shown on the ControlCOCKPIT display (T2, T3, etc.):




In the AtmoCONTROL software, the temperature sensor values are represented as additional coloured lines, also designated as T2, T3, etc.:



## HPP constant climate chamber with compressed air dehumidification


In this special variation, the temperature-humidity working area of the appliance is extended by the use of external compressed air. This makes it possible to dehumidify the air in the appliance even at 10 °C down to 10 % relative humidity:



In order to do this, connect a suitable pressure hose at the rear of the chamber (see picture) to an external compressed air supply (max. 5 bar).

**Only oil-free compressed air may be used.**

If 10 % humidity is not attained at 10 °C, the compressed air is not, or not sufficiently, pre-dehumidified. In this case it is possible to dehumidify the compressed air before it is fed into the chamber via a maintenance unit connected upstream, available from MEMMERT. Several chambers can be connected to this maintenance unit.




## Gas flushing

### Description

When equipped with gas flushing, gas can be flushed through the appliance. The gas flows in through a ball valve on the upper right and is channelled out through a second ball valve on the bottom left. The ball valves have a 3/8" internal thread to connect them to the system.

At the inlet valve (1), standard gas bottles with pressure-relief valve can be connected (maximum connection pressure 0.5 bar). Open the outlet valve (2) before injecting gas. Ensure that there is no overpressure in the appliance. The released gas must be channelled out.



Appliance with gas flushing (schematic diagram)

- 1 gas inlet
- 2 gas outlet

## Safety regulations

Observe the following special precautionary measures and safety regulations for appliances with gas flushing:



### **Warning!** **Danger of explosion and poisoning!**

- Only inject non-combustible, non-flammable, non-explosive, non-toxic and non-corrosive gases into the appliance.
- Always close the pressure-relief valve at the gas bottle and ball valves if the appliance is not in operation.
- Do not leave the appliance door open while gas is flowing in.
- Always keep the outlet valve open while injecting gas.
- Do not operate the appliance without ventilation at the outlet valve.
- Read the safety notes and instructions of the gas supplier.



## Handling

### Operation

1. Put the appliance into operation.
2. Open the outlet valve (2) on the bottom left of the appliance.
3. Open the gas bottle (max. 0.5 bar).
4. Open the inlet valve (1).

### Ending operation

1. Close the gas bottle.
2. Close the inlet valve (1).
3. Close the outlet valve (2).
4. Switch off the appliance.
5. Ventilate the appliance (open the door).









**memmert**  
Experts in Thermostatics



D24124 | Date 11/2016

englisch

Memmert GmbH + Co. KG  
Willi-Memmert-Straße 90-96 | D-91186 Büchenbach  
Tel. +49 9122 925-0 | Fax +49 9122 14585  
E-Mail: [sales@memmert.com](mailto:sales@memmert.com)  
[facebook.com/memmert.family](https://facebook.com/memmert.family)  
Die Experten-Plattform: [www.atmosafe.net](http://www.atmosafe.net)