# AtmoSAFE

## Was gibt es Neues? *What is the news?*

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### Requested by our customers: Intermediate size UIS 550

*There are many new things at our company to announce to you. In the first edition of our new customer newspaper, AtmoSAFE, you will learn on pages 3 and 4 how our CO<sub>2</sub> incubator INCO2 is used in tests of drinking water, and what it achieves here for the safety of employees and samples. We inform you of changes in equipment fittings on page 2, incidentally revealing a few details of our current development projects.*

*The title page has been given over to a new appliance: our UIS model range – universal ovens, incubators and sterilisers – has a new sibling. We have listened to your numerous suggestions, and with the 550 range we are addressing the need for an appliance size lying in-between our most popular size, the 500, with a 108 litre chamber capacity and the two-door 600, with a 256 litre capacity. From mid-October the new intermediate size UIS 550, with a chamber capacity of 153 litres, can be ordered in the controller classes EXCELLENT and PERFECT. Delivery will start from January 2009.*

### Our appliances are 100% AtmoSAFE

*The new and further development of our temperature-control ovens are guided by a crystal-clear vision: Controlled atmospheres. Working to achieve the highest possible degree of precision, reliability and flexibility is our passion, and this will in future define product development in our company. We have found a word for the outcome of this work: AtmoSAFE.*

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### New structure in Sales and Marketing

*One of our objectives for the coming years: to get even closer to the customer, provide even better care with respect to sales, marketing and customer service. For this reason from 1st October 2008, we have consolidated the internal organisation units that are particularly close to our customers, in the new department of Sales and Marketing. Our technical service, together with our customer service, will in future be led jointly by Gabriele Hartmann and Heinz Bayer, Angelika Henneberg officially takes over the marketing department. Gabriele Hartmann was until now responsible for sales and distribution in French and Spanish speaking countries. Heinz Bayer has been head of the Technical Service department since 1996. Angelika Henneberg has looked after customer service for the markets in Asia and the Middle East and coordinated Memmert marketing activities.*

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### A genuine innovation awaits you! Waterbaths with Peltier cooling

*New developments are normally top secret for all companies. Today we are making an exception to this rule and let you have a glimpse behind the scenes in our development department. If everything runs according to plan, in January 2009 we will have completed development work on waterbaths with Peltier cooling technology as an option for the standard fitting. A genuine innovation with which we, as the only manufacturer on the market, can offer a cooling option for waterbaths without expensive and energy-intensive cooling units. This means that for the first time constant temperatures between 10 °C and 30 °C will also be possible, since until now the temperature had to be at least 5 °C above room temperature. Appliances of the most up-to-date generation already in use can be easily upgraded with snap-on-technology. The waterbath based on Peltier technology is the third appliance series from our company, after the IPP incubator and the HPP constant climate chamber, which works with this cost-saving technology. In view of rising energy costs and increased efforts to protect the environment, they point the way forward for temperature-control technology, not only from an economic and ecological perspective. A future that for us has long since begun.*

### Vacuum ovens now with aluminium thermoshelves

*From January 2009 there is also a change in the fittings for vacuum drying ovens. To achieve an even better temperature distribution on all separately heated levels in the chamber, the thermoshelves will in future be manufactured from anodised and surface-enhanced aluminium. For applications with highly-corrosive materials we still recommend that you choose stainless steel shelves, because of their even higher resistance to chemicals.*

### Steel grids of stainless steel will be standard in UIS ovens

*We are remaining loyal to high-grade stainless steel – even for the steel grids with which our temperature-control ovens will be equipped by default from January 2009. If you base your product development on the requirements of your customers, as we do, then you know that you will never reach the point at which you are satisfied. And that's a good thing. Because only if we are continuously in dialogue with our customers can we get even better. The demand for steel grids for our appliances has increased in recent months, so we have decided to make this change to the standard equipment. You can still order stainless steel shelves, of course, at no extra charge.*


#### Your advantages at a glance:

- Better use made of the overall chamber, since an equal temperature distribution is guaranteed despite being heavily loaded
- High-quality, corrosion-proof stainless steel, no chrome-plated steel
- Rounded edges
- With no extra charge for standard scope of delivery
- Stainless steel shelves can be optionally chosen at no extra charge

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*Introduced in 2008 at the Analytica: HPP constant climate chamber for long-term tests in industry and research. High energy-saving effects through Peltier technology for heating and cooling.*

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#### Your advantages at a glance :

- Greater usable surface
- Optimised temperature distribution
- With no extra charge for standard scope of delivery
- Optional: thermoshelves made of stainless steel can be chosen at no extra charge

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## Growing bacteria under „cleanroom conditions“


Hygiene is the top priority when testing drinking water for Legionnaire's disease

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*For German hospitals, nursing homes and homes for the elderly, swimming baths and other public institutions, regular tests by certified laboratories of drinking water for Legionnaire's disease are prescribed. A part of the water samples taken in the Nuremberg conurbation arrive at the team led by Peter Daum in the Municipal Water Treatment and Environmental Analysis Laboratory Nuremberg, Germany.*

Memmert CO<sub>2</sub> Brutschrank INCO2  
Memmert CO<sub>2</sub> Incubator INCO2

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### Increased measures against Legionnaire's disease

*In 1999 the Environmental Analysis Laboratory Nuremberg put into operation the first Memmert CO<sub>2</sub> incubator for growing legionellae. Almost 10 years later it is being joined by the latest INCO2 generation of appliances – and with good reason. With the increased awareness of the danger of these bacteria in warmed up water, with which humans come into contact through drinking, bathing, showering or for medical applications, the number of samples taken and tested has also increased. Coming into contact with the bacteria is not in itself dangerous; it is rather the breathing in of droplets of water, or aerosols, containing the bacteria, into the lungs that can lead to the life-threatening Legionnaire's disease. Legionellae prefer water temperatures between 25 °C and 45 °C and prosper in quiet places where there is no movement and turbulence. Numerous types can be found everywhere in natural lakes and rivers and get into the drinking water supply through ground water, but it is only our rising standard of living that has turned Legionella pneumophila in particular, which is responsible for an estimated 90% of all cases of Legionnaire's disease, into a worldwide hazard. Because where hygiene standards are poor, or where there are structural flaws in technical systems for the supply of warm water, it can find ideal conditions to propagate. Legionnaire's disease, a form of pneumonia, is therefore a genuine disease of affluence, the name of which goes back to an epidemic in 1976 in which 182 former American soldiers were taken ill, of whom 29 died.*

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#### Sterilisation protects employees and samples

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*Benedikt Schaefer from the German Federal Environmental Office in Bad Elster, a proven expert on drinking water and a member of the Water Standards Committee, points out the need for extremely careful hygiene measures during the legionellae tests. Two reasons are crucial for this recommendation: The health of employees should not be exposed to even the slightest risk of contaminated aerosols, at the same time the sample quality is ensured over the long incubation period of up to 10 days. Help in the form of being absolutely germ-free can only be provided in this situation by sterilisation. The chamber of the INCO2, including the ventilation system, the water trays and all the sensors, can be sterilised in a 4-hour programme at 160 °C. The infra-red sensor for the CO<sub>2</sub> measurement was designed by the Memmert engineers specially so that it can withstand these high temperatures without problem, and is thus sterilised.*

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


About one thousand water samples reach the Environmental Analysis Laboratory Nuremberg every year. An estimated 25% of these contain legionellae, whereby the DVGW spreadsheet W551, in which the procedure for taking samples is regulated, only suggests, outside high-risk areas in hospitals, that the test intervals should be more frequent if the legionellae concentration is more than 100 CFU (colony-forming units) per 100ml, and stipulates further measures and tests only when the concentration is more than 1000 CFU. For 7 to 10 days the samples are incubated at 36 °C ( $\pm 2$  °C), 2.5 % CO<sub>2</sub> and 95% relative humidity. Once samples have been taken, they can no longer be reproduced, of course, and for this reason the safety and reliability of the appliances is the highest priority. In the nine years in which the CO<sub>2</sub> incubator has been running for Peter Daum and his team, almost non-stop, there has not been a single breakdown, and the test of temperature precision, performed twice a year by the internal quality assurance team, always yielded optimal values.

All Memmert appliances have an acoustic alarm, and this clearly audible warning signal is set off if the door is opened for too long or if there is a failure of the central supply, such as CO<sub>2</sub> gas, power supply, etc. One of the numerous safety functions allowing the employees in Nuremberg to quietly concentrate on their work. Just like the clear presentation of the current states for humidity, temperature and CO<sub>2</sub> in the display, which continuously shows the correct progression of the incubation in the laboratory – with the newest generation of CO<sub>2</sub> incubators from the company Memmert, down to the level of the water container. To ensure that the samples do not dry out, the electronic control of the INCO2 regulates not only temperature and CO<sub>2</sub> content, but also the relative humidity. Ideally the humidity content in the chamber should always be between 90% and 95% relative humidity.

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### Electropolishing: Small cause, large effect

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Progressive flattening of a metal surface through electropolishing

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The surface roughness of conventional stainless steel can neither be felt with the fingers nor seen with the bare eye. And yet in some circumstances they can provide a large surface for germs only a few micrometers in length. For working with biological materials the stainless steel chambers of all MEMMERT CO<sub>2</sub> incubators INCO 2 are therefore additionally smoothed by electrolysis for working with biological materials – a procedure that is standard in the manufacture of medical instruments, in order to remove the micro roughness on the metallic surface. At the same time the corrosion resistance of the metal is increased by this procedure. For the employees in the laboratory, cleaning and disinfection are also simpler due to the extremely smooth surfaces, and considerably more effective with respect to hygiene.

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### Humidity chamber HCP as an alternative to the INCO2

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Even if the currently valid DIN EN ISO 11731-2:2007 no longer stipulates that samples must be fumigated with CO<sub>2</sub>, it does announce a recommendation to the extent that carbon dioxide can promote the growth of certain legionellae. For laboratories that do not use CO<sub>2</sub>, the Memmert humidity chamber HCP represents a comparable alternative to the INCO2. This also offers all the advantages of the INCO2 mentioned above, such as sterilisability, safety and controllability - the rH concentration can be set via the active humidity control between 20% and 95%.

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## Future-compliant documentation

Asked about future developments in the field of standards and regulations for legionellae tests in the laboratory, Benedikt Schaefer mentions the subject of documentation. The validation of processes will become ever more significant, which is why suitable software solutions will have to be available to meet these requirements. Today the INCO2 is already fitted in its standard form with a whole series of convenient and time-saving functions and interfaces for programming, data storage and the documentation of processes.

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Memmert GmbH + Co. KG