



Fraunhofer

TESTED[®] DEVICE

Memmert GmbH + Co. KG
 Universal Oven UF1060
Report No. ME 1811-1085

Single product
Particle Emission

Qualification Certificate

This is to certify that the product mentioned above, provided by

Memmert GmbH + Co. KG
 Schwabach, Germany

has been awarded a Fraunhofer certificate TESTED DEVICE bearing the report number ME 1811-1085.

The Universal Oven UF1060 (color: stainless steel) was assessed in compliance with ISO 14644-14. When operated under the specified test conditions, it is suitable for use in cleanrooms fulfilling the specifications of the following Air Cleanliness Classes according to ISO 14644-1:

Test parameter(s)	Air Cleanliness Class
Outside: T = 300°C; flap = open; fan = 100 %	5
Inside: T = 300°C; flap = open; fan = 100 %	5
Inside: T = 300°C; flap = open; fan = 0 %	6
Overall result	6

DUPLICATE

This document only applies to the named product in its original state and is valid for a period of 5 years from the date the first document was issued. The document can be verified under www.tested-device.com.

Detailed information and parameters of the test environment can be found in the Fraunhofer IPA test report.

ME 1811-1085 Stuttgart, September 1, 2019
 Report No. first document Place, date of first document issued

-- --
 Report No. current document Place, current date

on behalf of 
 Dr.-Ing. Frank Bürger, Project Manager Fraunhofer IPA