



# Fraunhofer

## TESTED<sup>®</sup> DEVICE

Memmert GmbH + Co. KG  
Universal Oven UF1060  
**Report No. ME 1811-1085**

Single product  
**Particle Emission**

### Qualification Certificate

This is to certify that the product mentioned above, provided by

**Memmert GmbH + Co. KG**  
Schwabach, Germany

has been awarded a Fraunhofer certificate TESTED DEVICE bearing the report number ME 1811-1085.

The Universal Oven UF1060 (color: stainless steel) was assessed in compliance with ISO 14644-14. When operated under the specified test conditions, it is suitable for use in cleanrooms fulfilling the specifications of the following Air Cleanliness Classes according to ISO 14644-1:

Test parameter(s)	Air Cleanliness Class
Outside: T = 300°C; flap = open; fan = 100 %	5
Inside: T = 300°C; flap = open; fan = 100 %	5
Inside: T = 300°C; flap = open; fan = 0 %	6
<b>Overall result</b>	<b>6</b>

This document only applies to the named product in its original state and is valid for a period of 5 years from the date the first document was issued. The document can be verified under [www.tested-device.com](http://www.tested-device.com).

Detailed information and parameters of the test environment can be found in the Fraunhofer IPA test report.

DUPLICATE

ME 1811-1085  
Report No. first document

Stuttgart, September 1, 2019  
Place, date of first document issued

--  
Report No. current document

--  
Place, current date

on behalf of   
Dr.-Ing. Frank Bürger, Project Manager Fraunhofer IPA